



IN REPLY REFER TO:

# United States Department of the Interior

U.S. GEOLOGICAL SURVEY

Box 25046 M.S. 407

Denver Federal Center

Denver, Colorado 80225

## **NATIONAL WATER QUALITY LABORATORY TECHNICAL MEMORANDUM 1997.07**

March 26, 1997

To: Chief, Office of Water Quality  
Assistant Chief, Office of Water Quality  
Assistant Chief Hydrologist for Technical Support  
Regional Hydrologists  
Chief, NAWQA  
Chief, Office of Ground Water  
Assistant Chief, Office of Ground Water  
Chief, National Water Information System  
Chief, Office of Hydrologic Research  
Chiefs, Branches of Regional Research  
Area Hydrologists  
District Chiefs  
Regional Water Quality Specialists  
Assistant Regional Hydrologists for NAWQA  
District Water Quality Specialists  
Chiefs, NAWQA Study Units  
Chief, Quality Water Service Unit, Ocala  
Chief, Yucca Mountain Project  
QA Manager, Yucca Mountain Project  
Chief, Branch of Technical Development & Quality Systems  
Employees, National Water Quality Laboratory

From: Peter F. Rogerson, Chief  
National Water Quality Laboratory  
Branch of Analytical Services

Subject: Replacement Chemical Oxygen Demand Instrumentation Operational

Author: Glenda Brown, Inorganic Program (303) 467-8122 (gebrown@usgs.gov)

Revision: None

Supplemental: None

### **SCOPE**

Effective December 9, 1996, new instrumentation for Chemical Oxygen Demand (COD), Method I-3561-85, has been put into routine operation by National Water Quality Laboratory (NWQL)

personnel. Chemical methodology and reagents are unchanged from TWRI Method I-3561-85. To document the instrumentation changes, a new Lab Code (2144) and a new method code (340C) have been assigned. The method reporting limit remains at 10 mg/L. The bottle type and preservation remain the same as Lab Code 76 (LC0076).

## BACKGROUND

TWRI Method I-3561-85 for COD determinations has used the same instrumentation and technology for the past 30 years. This instrumentation has been replaced with a Hach COD digester and colorimeter system. The new system is considerably safer to operate than the old system and generates much less hazardous waste. The Hach COD system uses EPA-approved technology.

Because of instrumentation problems, any samples received by the laboratory after October 23, 1996, have been analyzed using the Hach COD system. Reruns will be analyzed by the Hach COD system. Samples received by the laboratory between October 23 and December 8, 1996, were held until the Hach COD system was approved for use.

## NEW INSTRUCTIONS AND WHAT TO EXPECT

Lab Code 76 is no longer valid; use Lab Code 2144 instead. Please discontinue use of Lab Code 76 and instead label the COD sample bottles with the new 2144 lab code. Requests for Lab Code 76 will be transferred to the new lab code. Continue to collect samples for COD analysis in a 125-mL glass bottle, preserve them with sulfuric acid to pH <2, and chill the samples.

Data quality achieved with the Hach COD system is equivalent to or better than that achieved with the retired COD system. Precision measured for the Hach COD system, expressed in relative standard deviation (RSD), is shown in the table below. RSD charts are available on request from the author at the NWQL.

The Hach COD system will provide NWQL customers with equivalent or better data in a more timely and safe manner.

Samples with concentrations between 1 and 10 mg/L COD may be analyzed on a custom basis. Contact Glenda Brown (gebrown@usgs.gov) for further details.

---

---

### PRECISION OF HACH COD SYSTEM

CONCENTRATION (mg/L)	PERCENT RELATIVE STANDARD DEVIATION
10-20	<20
20-50	<10
50-100	<5
>100	<3

---

---

Keywords: COD, chemical oxygen demand

Supersedes: None

Effect on Database: None

Distribution: See above plus the Netnews usgs.labnews & .water.quality; WRD Secretaries; Field and Project Offices; Hydrologic Technicians; and <http://wwwnwql.cr.usgs.gov/>